Dual Top Anchor System

■ Dual Top Anchor System
■ Sterile Dual Top Anchor System

INSTRUCTIONS FOR USE

Manufacturer:
JEILMEDICAL
702-703-704-705-706-804-805-807-812-815-ho, 55, Digital-ro 34-gil,
Guro-gu, Seoul, 08378, Rep. of Korea Tel 82 2 850 3500 Fax 82 2 850 3535
작성연월(Issue Date) : 2023.09 **Made in Korea**
Print Date : #027-41001-01

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English

[DESCRIPTION]

The Dual Top Anchor System was designed to provide a firm anchorage in the oral cavity for an orthodontic treatment. Orthodontic screws are made of a titanium alloy (ASTM F 136). The head of the Dual Top Anchor System screw is designed to apply to various orthodontic tools. The sizes of the Dual Top Anchor System are diverse enough to satisfy various clinical cases. This system includes screw blocks and tools such as drill bit, driver shaft and hand body.

[DIRECTIONS FOR USE]

1. Screw insertion without drilling

- ① Select a screw and a driver.
- ② Administer a topical anesthetic to a graft site.
- ③ Attach the selected screw to the driver.
- ④ Place the screw straight on the mucosa, and then turn the handle as pressing it so that the tip of the screw can dig a hole in the surface cortical bone. (While inserting the screw, turn it slowly to prevent the driver axis from shaking.)

- ⑤ Stop inserting when the head at the lower end of the screw meets the mucosa, and remove the driver from the screw.

2. Screw insertion after drilling

- ①,② The same as above.
- ③ Attach pilot drill to a hand piece, and start drilling. (1,300 r.p.m or less. The drilling area should be cooled by physiologic saline.)
- ④ Attach the selected screw to the driver.
- ⑤ Turn the handle, fitting the hole in the cortical bone. (While inserting the screw, turn it slowly to prevent the driver axis from shaking.)
- ⑥ Stop inserting when the head at the lower end of the screw meets the mucosa, and remove the driver from the screw.

3. Screw removal

- ① Select a driver that fits the screw.
- ② Set the selected driver on the screw head, and then remove the screw by turning it slowly in the opposite direction of the screw insertion.

[INDICATIONS FOR USE]

Dual Top Anchor System is used to provide a fixed anchorage point for attachment of orthodontic appliances to facilitate the orthodontic movement of teeth. It is used temporarily and is removed after orthodontic treatment has been completed. Screws are intended for single use only. For use on adults over the age of 12.

[**MATERIAL**] Titanium Alloy (ASTM F 136)

[WARNING]

- U.S. Federal Law restricts this device to sale by or on the order of a licensed surgeon or other qualified medical professional.

Selection of screw size must be carefully considered by the operating surgeon and should take into consideration the quality of bone, bone type, functional loads exerted on bone(s), and post-operative patient compliance.

- For sale and effective use of [Dual Top Anchor System], the surgeon should have specific training, experience, and thorough familiarity with the use of rigid fixation and techniques.
- The surgeon must exercise reasonable judgment when deciding which screw type to use for specific indications.
- Dual Top Anchor System are not intended to endure excessive abnormal functional stresses.
- Use of an undersized screw in areas of high functional stress may lead to implant fracture and failure.

[CAUTION]

- Fixing the product incompletely may cause damage or breakage of the product, and the patient may feel a foreign body sensation.
- Always follow appropriate safety precautions.
- Select the appropriately sized screw for the patient.
- Responsibility for proper selection of patients, adequate training, experience in the choice and placement of plate & screw and the decision to leave or remove plate and screw postoperatively, rests with the surgeon.
- Once applied, never reuse this device.
- When engaging the screw, axial pressure of the driver

- shaft into the screw head must be adequately applied to ensure that the blade is fully inserted into the screw head. This results in proper axial alignment and full contact between driver shaft and screw, minimizing the risk of round-out. Otherwise there will be an increased risk of mechanical damage to the implant or the driver shaft.
- Pilot drilling is required, if the bone density of the patient is too high or too thick.
- Careful handling and storage of the product is required. Store in a dry place at room temperature.

[PRECAUTIONS]

- The surgeon should inform the patient of risks associated with surgery.
- Inspect each device to ensure they are not bent or damaged.
- Before use, check that the specifications of the product matches those on the label.

[POSSIBLE ADVERSE EFFECTS]

- In many cases, adverse results may be clinically related rather than implant related.
- Osteoporosis, inhibited revascularization and poor bone formation can cause loosening, bending, cracking and fracturing of this device or premature loss of fixation with the bone, leading to nonunion.
- Inadequate arrangements may cause procedure failure.
- The product should not be used on patients sensitive to metals such as titanium, as it may cause boils and inflammatory reactions.

[CONTRAINDICATION]

- Not for use in cases of active or suspected infection or in patients previously sensitized to titanium.
- Not for use in patients exhibiting disorders which would cause the patient to ignore the limitations of rigid fixation plate and screw implants.

[CLEANING]

Dual Top Anchor System's accessories (surgical instruments) can be reused after cleaning and sterilizing. They must be cleaned before reuse. Following instructions are recommended for cleaning;

- 1) Pre-cleaning
 - ① Disassemble the device where possible.
 - ② Remove gross soil using paper wipes and solution of cleaning agent.
 - ③ Immerse the device in solution of cleaning agent for the time recommended by the detergent's manufacturer.
 - ④ Using suitable brushes (never metal brushes or steel wool) cleaning the device thoroughly.
 - ⑤ Rinse in running water until all traces of cleaning solution are removed.
 - ⑥ Visually inspect for any remaining soil and repeat the steps above if necessary.

2) Cleaning

- ① Immerse the device completely in solution of cleaning agent and activate the ultrasonic bath at the concentration and temperature specified in the detergent manufacturer's instructions.
- ② Follow the same as the clause ④ ~ ⑥ of instructions for pre-cleaning.

[STERILIZATION]

• For Non-sterile

Dual Top Anchor System must be sterilized before use. Following instructions are recommended;

- ① Place the devices in the appropriate block (kit or tray) using forceps and/or powder-free gloves to avoid contamination and any other negative effect on the surface of device.
- ② Wrap the block with a surgical drape.
- ③ Sterilize in the autoclave validated and maintained in accordance with ISO 17665 and ANSI AAMI ST79. Following parameters are validated in accordance with ISO 17665-1 and recommended for sterilization;

Cycle	Temperature	Exposure Time	Load Characteristics	Drying Time
Gravity	132°C	15 min.	Wrapped*	30 min.
Pre-vacuum	132°C	4 min.	Wrapped*	30 min.

- * In the case of load characteristics, we recommend usage of an FDA cleared wrap to ensure that the device is actually sterile prior to implantation.
- ④ Caution: Screws are disposable (single use only). Do not reuse. Jeil Medical Corporation is not responsible for all problems caused by reuse.

• For sterile

Sterile Dual Top Anchor System is supplied by gamma sterilization. Expiration date : 5 years from Manufacture date(Refer to Label)

[PACKAGING]

Dual Top Anchor System devices are packaged individually or set. This device manufactured and sold by Jeil Medical Corporation.

[SYMBOL DESCRIPTIONS]

 REF	Catalogue number	 LOT	Batch code
 Do not re-use		 Do not re-sterilize	
 Date of manufacture		 Manufacturer	
 Non-Sterile		 Caution	
 Use-by Date		 Sterilized using irradiation	
 Consult instruction for use			
 Do not use if package is damaged			