JEILMEDICAL

CE
1639

SMARTO

Instruction For Use

Issue Date: 2023. 01

Print Date:

#120-41004-01

Table of Content

English	03
Deutsch	12
Español	22
Ελληνικά	31
Italiano	40
Magyar nyelv	49
Türkçe	58

[DESCRIPTION]

The SMARTO is a sterile battery powered screwdriver. The device includes a DC motor, battery, switch and holder for the rotation attachment. It is disposable.

[SPECIFICATION]


Model name (111-ED-)	030	031	050	051	052
Torque (Ncm \pm 10%)	24	30	45	40	35
RPM (\pm 10%)	160	250	210	165	110

* Type of protection against electric shock:

- Internally powered equipment
- Degree of protection against electric shock
: Type B applied part

Holder Cap(Type B applied part)

Driver Shaft(Type B applied part)

- This equipment is not suitable for use in the presence of flammable anesthetics or oxygen

[INDICATION FOR USE]

The SMARTO is used to drive screws and drilling in conjunction with Orthopedic, Craniomaxillofacial, Neuro surgical procedures. It is supplied sterile and is for single use only. Recommended use for patients who are 18 and older.

[DIRECTIONS FOR USE]**1. Point to Check Before Using**



- ① Inspect package and device for signs of damage. Do not use if package is open or torn.
- ② This device is supplied sterile.

2. Operating Condition


- ① Temperature: 10~35 °C
- ② Humidity: 20~85% RH
- ③ Atmospheric: 500~1,060hPa
- ④ Altitude: 0~2,000m
- ⑤ Illuminance: 100~500lx

3. How To Use


- ① Insert the sterilized driver shaft the one most suitable for the surgery. (SMARTO is to be used with JEIL accessories only.)


- ② Set the motor in motion by pushing the button. (CLOCKWISE / COUNTERCLOCKWISE)

- ③ Set the Lock switch to use in manual mode (LOCK), when you need more torque than driver torque.
(Applied model: 111-ED-031, 050, 051, 052)

[DC-Powered]

[Manual]

- ④ After operation, dispose of waste materials as biohazard or according to local medical waste disposal requirements.

[PRECAUTIONS]

- ① To avoid any accidents when using this product, the user must be fully attentive for the safety of both the user and the patient.
- ② Only qualified medical practitioners with the necessary medical and technical licenses are allowed to use the product.
- ③ Check the packaging condition for signs of damage.
- ④ Before starting any treatment on a patient, test this product first by rotating it to check the product is operating smoothly.
- ⑤ If any abnormality is detected with the product before use (such as vibration, noise, overheating, rotation problems etc.), immediately contact the distributor. Replace with a new Smarto.
- ⑥ The driver shaft must be cleaned and sterilized before use. (A dirty driver shaft could infect patients with germs. A motor stall could also occur due to dust or impurities.) SMARTO is to be used with JEIL accessories only.
- ⑦ The product must be off while the driver shaft is being replaced. An accident could occur if the driver shaft is replaced while the driver is rotating.
- ⑧ Do not attempt to disassemble or modify the product. Only qualified engineers of the manufacturer can repair or modify the product.

[CONTRAINDICATION]

Do not recommend for use in patients with low bone density.

[WARNING]**1. SAFETY INSTRUCTION**

- ① Do not operate both switches (CLOCKWISE/COUNTERCLOCKWISE) at the same time.
- ② If any abnormality is detected in the product during use (such as vibration, noise, overheating, rotation problems etc.), stop using the product immediately and contact the distributor.
- ③ Products intended for single use only must not be re-used. Devices labeled for single-use only should never be reused. Re-use of these devices may cause harmful effects to patients. Re-use or reprocessing (e.g. cleaning and resterilization) may compromise significant degradation in device performance, crossinfection, and contamination. JEILMEDICAL Corporation does not take any responsibilities for any problems caused by reuse.
- ④ Be aware of overheating when continuously used.

2. DISPOSAL

- ① After operation, dispose of waste materials according to disposal method of medical device.
- ② Do not disassemble or modify the product at your own discretion. Only qualified engineers of the manufacturer can repair or modify the product.

3. STORAGE

- ① Keep the product away from direct sunlight, high temperature, flammable liquids or gases.
- ② Temperature: 10~35°C
- ③ Humidity: 20~85% RH
- ④ Atmospheric: 500~1,060hPa

[MATERIAL] ABS etc.

[HANDLING & STERILIZATION] This device is supplied in sterile condition.

[CAUTION] This device is restricted to sale, distribution and use by or on the order of a medical practitioner licensed to do so, and use by other person is prohibited.

[PRECAUTIONS FOR ELECTROMAGNETIC WAVES]

This equipment has been tested and found to comply with the limits for medical devices in IEC/EN 60601-1-2:2014. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation.

This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving device.
- Increase the separation between the equipment e equipment.
- Connect the equipment into an outlet on a circuit different from that to which the other device(s) are connected.
- Consult the manufacturer or field service technician for help.

[MANUFACTURE'S DECLARATION]

Manufacturer's declaration - electromagnetic emission

JEILMEDICAL, Sterile Single Use Battery Powered Surgical Handpiece are intended for use in the electromagnetic environment specied below. The customer or the user of this should assure that it is used in such an environment.

Emission test	N/A	N/A
RF emissions CISPR 11	Group 1	The Model 111-ED-050 RF energy only for its internal function. Therefore. Its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
	Class B	
Harmonics emission IEC 61000-3-2	N/A	N/A
Voltage fuctuation IEC 61000-3-3		

Manufacturer's declaration - electromagnetic immunity

JEILMEDICAL, Sterile Single Use Battery Powered Surgical Handpiece are intended for use in the electromagnetic environment specified below. The customer or the user of this should assure that it is used in such an environment.

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic Environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-2	6 kV Contact 8 kV Air	6 kV Contact 8 kV Air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast Transient/burst IEC 61000-4-4	N/A	N/A	N/A
Surge IEC 61000-4-5	N/A	N/A	N/A
Power frequency (50/60Hz) Magnetic eld IEC 61000-4-8	3.0 A/m	3.0 A/m	Power frequency magnetic elds should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Voltage dips, short Interruptions and Voltage variations on power supply input lines IEC 61000-4-11	N/A	N/A	N/A
Conducted RF IEC 61000-4-6	N/A	N/A	N/A

Radiated RF IEC 61000-4-3	3 V/m 80.0 MHz to 2.5 GHz	3 V/m 80.0 MHz to 2.5 GHz	Field strengths outsourced RF transmitters, as determined by an electromagnetic site survey, should be less than 3V/m.a Interference may occur in the vicinity of equipment marked with the following symbol: 
------------------------------	---------------------------------	---------------------------------	--

Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and the JEILMEDICAL, Sterile Single Use Battery Powered Surgical Handpiece.

The Model 111-ED-050 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The user of the 111-ED-050 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the 111-ED-050 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output Power (W) of transmitter	Separation distance (m) according to frequency of transmitter		
	150 kHz bis 80 MHz $d = \left[\frac{35}{\sqrt{f_1}} \right] \sqrt{P}$	80 MHz bis 800Hz $d = \left[\frac{35}{\sqrt{f_1}} \right] \sqrt{P}$	800 MHz bis 2.5 GHz $d = \left[\frac{7}{\sqrt{f_1}} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

[EXPIRY DATE] Up to three years from the date of manufacture.

[Sterile by gamma irradiation]

[SYMBOL DESCRIPTIONS]

<i>Symbol</i>	<i>Description</i>	<i>Symbol</i>	<i>Description</i>
	Catalogue number		Batch code
	Do not re-use		Do not re-sterilize
	Date of manufacture		Manufacturer
	Use-by Date		Sterilized using irradiation
	Caution		Consult instruction for use
	Do not use if package is damaged		Type B applied part
	Clockwise		Counterclockwise
UNLOCK	DC-Powered mode	LOCK	Manual mode
	Authorized representative in the European Community		EMERGO EUROPE Westervoortsedijk 60, 6827 AT Arnhem, The Netherlands

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