ARIX System Bone Screw / Bone Plate

- ARIX Hand System
- ARIX Foot System
- ARIX Wrist System
- ARIX Ankle System
- ARIX Ankle Fibula Hook Plate System
- ARIX Ankle Distal Tibia System
- ARIX Diaphysis System
- ARIX Humerus System
- ARIX Clavicle System
- ARIX Cannulated Screw System
- ARIX Elbow System
- ARIX Small Fragment System

INSTRUCTIONS FOR USE

JEILMEDICAL Manufacturer: Jeil Medical Corporation
702-703-704-705-706-804-805-807-812-815-ho
55, Digital-ro 34-gil, Guro-gu, Seoul, 08378, Rep of Korea Tel 82 2 850 3500 Fax 82 2 850 3535

#467-41002-10 **Made in Korea**

작성일월 (Issue Date) : 2024.03

Print Date :

English

[DESCRIPTION]

The ARIX Hand / Foot / Wrist / Ankle / Ankle Fibula Hook Plate / Ankle Distal Tibia / Diaphysis / Humerus / Clavicle / Cannulated Screw / Elbow / Small Fragment System is rigid fixation consisting of plates and screws in various configurations, shapes and sizes, as follows:

System	Screw	Plate
ARIX Hand System	12/15/18-HC Series	H1 Series
	20/23/25-HC Series	H2 Series
	15-HC Series 15L-HF Series 20/20L-HF Series 23-FC Series	H1L Series H2L Series
ARIX Foot System	23/30/45/60-CS Series 15/23/30-CH Series	N/A
	28-FC Series 28L-HF Series 28/28L-SO Series	F28 Series
	35-FC Series 35L-HF Series 35/35L-SO Series	F35 Series Calcaneal Series
	28/28L-SO-L Series 35/35L-SO-L Series	F35 Series Calcaneal Series
	23/23L-SO-S Series 28/28L-SO-S Series 20-AO Series 27-AO Series	Mid-Fore Series Wedge Series
ARIX Wrist System	161.0325 Series 25L/25-HF Series 25C/25P-HF Series 25/25L-SO Series 28/28L-SO Series 28C-HF Series 28P-HF Series 28/28L-HF Series 28V-HF Series 20/20L-SO Series 15L-HF Series 25P-SO Series	25-DVRA Series 25-DLUL Series 25V-DLUL Series 25-DDRA Series 20-DUHK Series 25-CAFU Series 25R-DVRA Series 25J-DVLA Series
	25/25L-SO Series 28/28L-SO-S Series 35-FC Series 35L-HF Series 40-CS Series 35-SO-L Series 35L-SO-L Series	35-DLFI Series 35-DLST Series 35V-DLF2 Series 35V-DLF3 Series
ARIX Ankle System	35-SO-L Series 35L-SO-L Series 35-FC Series 35L-HF Series	35-DLFH Series 35V-DLFH Series
ARIX Ankle Distal Tibia System	28-SO-L Series 28L-SO-L Series 35-SO-L Series 35L-SO-L Series 40-SA Series 40-SA Series	35-DLTB Series 35-DMTB Series 35-DATB Series
	50-SO Series 50L-SO Series 55-SA Series 55L-SA Series	50-AAFU Series
ARIX Diaphysis System	35-SO-L Series 35L-SO-L Series	35-DIST Series 35V-DIST Series
ARIX Humerus System	35-SO-L Series 35L-SO-L Series 40L-SA Series 45L-CO Series	35-PLHU Series
ARIX Clavicle System	35-SO-L Series 35L-SO-L Series 25L-SO Series	35-SMCL Series 35-SLCL Series 35-HPCL Series
ARIX Cannulated Screw System	23/30/45/60-CS Series 23/30-CH Series 40-CO Series 65-CO Series 73-CO Series	N/A
ARIX Elbow System	20/20L-SO Series 25/25L-SO Series 28/28L-SO Series 35-SO-L Series 35L-SO-L Series	20-PLRA Series 25-PLUL Series 35-PLUL Series 35-DMHU Series 35-DLHU Series 35-DPHU Series
	28/28L-SO Series	25-SFST Series 25-SFTR Series 25-SFTS Series

The ARIX Hand / Foot / Wrist / Ankle / Ankle Fibula Hook Plate / Ankle Distal Tibia / Diaphysis / Humerus / Clavicle / Cannulated Screw / Elbow / Small Fragment System is made of Pure Titanium and Titanium Alloy (Ti-6AL-4V), which meet ASTM F67, Standard Specification for Unalloyed Titanium for Surgical Implant Applications, and ASTM F136, Standard Specification for Wrought Titanium-6 Aluminum-4 Vanadium ELI (Extra Low Interstitial) Alloy for Surgical Implant Applications, which are widely used for surgical implants with well-known biocompatibility.

[DIRECTIONS FOR USE]

The ARIX Hand / Foot / Wrist / Ankle / Ankle Fibula Hook Plate / Ankle Distal Tibia / Diaphysis / Humerus / Clavicle / Cannulated Screw / Elbow / Small Fragment System Plate & Screw offers the operating surgeon unique features and user friendly capabilities. A step by step explanation of how to use the system follows:

1.Selecting Plate Configurations

The choice of which plate or screw size to use is at the discretion of the operating surgeon.

2.Plate Handling

The ARIX Hand / Foot / Wrist / Ankle / Ankle Fibula Hook Plate / Ankle Distal Tibia / Diaphysis / Humerus / Clavicle / Cannulated Screw / Elbow / Small Fragment System Plate & Screw is grasped, using the forceps supplied with the ARIX Hand / Foot / Wrist / Ankle / Ankle Fibula Hook Plate / Ankle Distal Tibia / Diaphysis / Humerus / Clavicle / Cannulated Screw / Elbow / Small Fragment System.

3.Plate Bending

The ARIX Hand / Foot / Wrist / Ankle / Ankle Fibula Hook Plate / Ankle Distal Tibia / Diaphysis / Humerus / Clavicle / Cannulated Screw / Elbow / Small Fragment System Plate can be bent and contoured using a plate bender. Exercised care to avoid distorting the screw hole.

4.Drilling Pilot Holes

The ARIX Hand / Foot / Wrist / Ankle / Ankle Fibula Hook Plate / Ankle Distal Tibia / Diaphysis / Humerus / Clavicle / Cannulated Screw / Elbow / Small Fragment System's screws are self-tapping. The cortical tip design and sharp threads allow the screws to self-tap into bone. All self-tapping screws require that a pilot hole be made in the bone prior to screw insertion. Drilling a pilot hole reduces the chance of fracturing bone and reduces the level of torque required to drive screws. To reduce the potential for thermal necrosis, all pilot drills should be used at minimal RPM. Proper pilot drilling technique indicates that the drill point do the cutting. Fault drilling causes that screws will not fit tightly into the bone or even strip in the bone. The surgeon must use drill fit the diameter of screws.

5.Plate and Screw fixation

To fixate plate and screw, the surgeon uses driver shafts. Insert the driver shaft in the handbody. Use it to pick up a screw then rotate it to fixate the screw in the surgical site. A driver shaft features a cruciform tip that friction fits into either the screw heads. A driver shafts or other instruments are separated from the screwdriver body by pulling back on the collar at the distal end of the screwdriver body.

*Note: Smooth peg screw is intended to support subchondral bone as a buttress rather than fixation of bone fracture. Depending on the type of bone fracture, Smooth Peg Screw can be used optionally at the discretion of the operating surgeon.

6.Plate & Screw Remove

In instances where ARIX Hand / Foot / Wrist / Ankle / Ankle Fibula Hook Plate / Ankle Distal Tibia / Diaphysis / Humerus / Clavicle / Cannulated Screw / Elbow / Small Fragment System removal is required, after normal surgical exposure the ARIX Hand / Foot / Wrist / Ankle / Ankle Fibula Hook Plate / Ankle Distal Tibia / Diaphysis / Humerus / Clavicle / Cannulated Screw / Elbow / Small Fragment System Plate & Screw may be dislodged by using a driver shaft and other instruments. If ARIX Hand / Foot / Wrist / Ankle / Ankle Fibula Hook Plate / Ankle Distal Tibia / Diaphysis / Humerus / Clavicle / Cannulated Screw / Elbow / Small Fragment System Plate & Screw removal are required, they should be removed, inventoried and discarded. Do not reuse. The surgical site is now re-sutured.

7.Cannulated Screw(Lag Screw Technique)(When applicable)

-Insert a guide pin through the drill guide to the appropriate depth. Check placement of the wire under image intensification. Remove the drill guide.

-In areas where soft tissue coverage is minimal, use the cannulated countersink to create a recess for the screw head. Countersinking will also facilitate screw insertion.

-Slide the tapered end of the cannulated screw measuring device over the guide pin, down to the bone. Read the scale at the end of the guide pin to determine appropriate screw length. This reading will place the screw 5mm short of the tip of the guide pin, do not remove the guide pin in the bone during screw insertion.

-Place the appropriate length screw over the guide pin. Use the cannulated driver shaft and handbody to insert the screw. Remove and discard the guide pin.

[INDICATIONS]

System	Intended for use
ARIX Hand System (K131566)	The ARIX Hand System is intended for use in internal fixation of the bones of hand and wrist. Examples of these procedures may include but are not limited to replantation, lag screw techniques, joint fusions, corrective osteotomies, and the treatment of fractures.
ARIX Foot System (K131311)	The ARIX Foot System is intended for use in internal fixation, reconstruction or arthrodesis of small bones including the fore, mid- and hind foot ankle. Examples of these procedures may include but are not limited to replantation, lag screw techniques, joint fusions, corrective osteotomies, and the treatment of fractures.
ARIX Wrist System (K151468, K191972)	The ARIX Wrist System(Radius) is intended for use in forearm fractures, osteotomies and arthrodesis. The ARIX Wrist System(Ulna) is intended for fractures and osteotomies, in particular for the ulna.
ARIX Ankle System (K152158, K231441)	The ARIX Ankle System (Fibula) is intended for use in internal fixation of the distal fibula. The ARIX Ankle System is intended for use in internal fixation of the distal fibula.
ARIX Ankle Fibula Hook Plate System (K170979)	The ARIX Ankle Fibula Hook Plate System is intended for use in internal fixation of the distal fibula.
ARIX Ankle Distal Tibia System (K170313, K231887)	The ARIX Ankle Distal Tibia System is intended for fixation of complex intra- and extra-articular fractures and osteotomies of the distal tibia including distal tibia fractures in combination with diaphyseal fracture.
ARIX Diaphysis System (K171285)	The ARIX Diaphysis System is intended for fixation of fractures, osteotomies, and non-unions of the clavicle, scapula, olecranon, humerus, radius, ulna, pelvis, distal tibia, fibula, particularly in osteopenic bone.
ARIX Humerus System (K172008)	The ARIX Humerus System is indicated for fractures, fracture dislocations, osteotomies and non-unions of the proximal humerus, particularly in osteopenic bone.
ARIX Clavicle System (K180972)	The ARIX Clavicle System is indicated for fixation of single, segmental and comminuted fractures, osteotomies, mal-unions, and non-unions of the clavicle.
ARIX Cannulated Screw System (K192417, K233912)	The ARIX Cannulated Screw System is intended for use in internal fixation, reconstruction, or arthrodesis of small bones including the fore, mid- and hind foot and ankle. Also, the ARIX Cannulated Screw System is intended for use in internal fixation of the bones of hand and wrist. Examples of these procedures may include but are not limited to replantation, lag screw techniques, joint fusions, corrective osteotomies, and the treatment of fractures. The ARIX Cannulated Screw System is indicated for fracture fixation, fusions, osteotomies, nonunions, and malunions of fragments of long bones of the pelvis, the femur, the tibia, the humerus, and the foot. The system is not intended for use in the spine.
ARIX Elbow System (K201656)	The ARIX Elbow System is intended for fractures and osteotomies of proximal radius, ulnar olecranon and distal humerus. The ARIX Proximal Radius Plate is intended for use in proximal radial fractures and osteotomies. The ARIX Olecranon Plate is intended for use in particular for ulna fractures and osteotomies. The ARIX Distal Humerus Plate is intended for use in distal humerus fractures, osteotomies and non-unions.
ARIX Small Fragment System (K221920)	The ARIX Small Fragment System is intended for fixation of fractures, osteotomies, nonunions, replantations, and fusions of small bones and small bone fragments. Examples include, but are not limited to, the hand, wrist, foot, and ankle.

[WARNING]

- U.S. Federal Law restricts this device to the sale by or on the order of a licensed physician or other qualified medical professional.
- Use of undersized plates and screws in areas of functional stress may lead to implant fracture and failure.

- Plates, Screws and other devices of dissimilar material should not be used together in or near the implant site. dissimilar metals in contact with each other can accelerate the corrosion process due to galvanic corrosion effects. Do not mix different metal implants in the same construct.
- Multiple bending may weaken the plat and could result in implant fracture and failure.
- Use of screws in high density bone may lead to implant fracture or failure upon insertion.
- The ARIX Hand / Foot / Wrist / Ankle / Ankle Fibula Hook Plate / Ankle Distal Tibia / Diaphysis / Humerus / Clavicle / Cannulated Screw / Elbow / Small Fragment System is recommended for use in patients with sufficient bone quality to sustain effectiveness and benefits of rigid fixation.
- It is recommended to remove any fractured implants from patient during surgery. If unable to remove, notify patient.
- Use of excessive torque during insertion of screws may lead to implant failure.

[MR SAFETY INFORMATION]

- The ARIX Hand / Foot / Wrist / Ankle / Ankle Fibula Hook Plate / Ankle Distal Tibia / Diaphysis / Humerus / Clavicle / Cannulated Screw / Elbow / Small Fragment System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of ARIX Hand / Foot / Wrist / Ankle / Ankle Fibula Hook Plate / Ankle Distal Tibia / Diaphysis / Humerus / Clavicle / Cannulated Screw / Elbow / Small Fragment System in MR environment is unknown. Scanning a patient who has this device may result in patient injury.

[CAUTION]

- Do not attempt a surgical procedure with faulty, damaged, or suspect ARIX Hand / Foot / Wrist / Ankle / Ankle Fibula Hook Plate / Ankle Distal Tibia / Diaphysis / Humerus / Clavicle / Cannulated Screw / Elbow / Small Fragment System. Inspect all components preoperatively to assure utility. Alternate fixation methods should be available intraoperatively.
- Fixing the product incompletely may cause damage or breakage of the product, and the patient may feel a foreign body sensation.
- Always follow appropriate safety precautions.
- Select the appropriately sized plate and screw for the patient.
- Responsibility for proper selection of patients, adequate training, experience in the choice and placement of plate & screw and the decision to leave or remove plate and screw postoperatively, rests with the surgeon.
- Once applied, never reuse this device.

- It is recommended that the implant be removed after use for up to 16 weeks, which is the period of completion of bone healing, unless medical attention is required. (However, the removal period may be extended if there is a special opinion of the surgeon.)
- Delayed healing, nonunion or subsequent bone resorption or trauma may cause excessive stress on this device and result in loosening or fracture.
- Careful handling and storage of the product is required. Store in a dry place at room temperature.
- **Disposal:** This Surgical Instrument is recommended to use up to 10 times, provided that Jeil medical-recommended conditions of use are followed. The users can use more devices based by their judgment. It provides a manual for checking the instrument condition, and the company is not responsible for using the instrument according to the user's judgment.

[PRECAUTIONS]

- The physician should inform the patient of risks associated with surgery.
- Inspect each device to ensure they are not bent or damaged.

[POSSIBLE ADVERSE EFFECTS]

- In most cases, potential complications have a clinical source as opposed to arising from the ARIX Hand / Foot / Wrist / Ankle / Ankle Fibula Hook Plate / Ankle Distal Tibia / Diaphysis / Humerus / Clavicle / Cannulated Screw / Elbow / Small Fragment System
- These include among other things;
- Loosening of the implant from insufficient fixation
 - Hypersensitivity to metal or allergic reactions
 - Bone necrosis, osteoporosis, insufficient revascularization, bone resorption and poor bone formation that can cause premature loss of fixation
 - Soft tissue irritation and/or nerve damage due to surgical trauma
 - Early or late infection, both superficial and deep
 - Elevated fibrotic tissue reaction around the surgical area
 - Complications in implant from improper explanation of the implant

[CONTRAINDICATIONS]

- Use of the ARIX Hand / Foot / Wrist / Ankle / Ankle Fibula Hook Plate / Ankle Distal Tibia / Diaphysis / Humerus / Clavicle / Cannulated Screw / Elbow / Small Fragment System is contraindicated in following cases;
- Active or suspected infection or in patients who are immunocompromised.
 - Patients previously sensitized to titanium.
 - Patients with certain metabolic diseases.
 - Patients exhibiting disorders which would cause the patient to ignore the physician pre and/or post-operative instruction and limitations of internal rigid fixation implants.

[NON-STERILE DEVICES-IMPLANT]

Products supplied in a non-sterile condition must be steam sterilized prior to surgical use. Prior to steam sterilization, place the product in an approved sterilization wrap or container.

[REUSABLE DEVICES - Instrument, Trays]

Reusable devices (Surgical instruments, Trays) can be reused after cleaning and sterilization.

[CLEANING & STERILIZATION]

- Cleaning
ARIX Hand / Foot / Wrist / Ankle / Ankle Fibula Hook Plate / Ankle Distal Tibia / Diaphysis / Humerus / Clavicle / Cannulated Screw / Elbow / Small Fragment System's accessories (surgical instruments) can be reused after cleaning and sterilizing.
They must be cleaned before reuse. Following instructions are recommended for cleaning;
※ Please Refer to the Cleaning and Sterilization Guide.

1) Pre-cleaning

- ① Disassemble the device where possible.
- ② Remove gross soil using paper wipes and solution of cleaning agent.
- ③ Immerse the device in solution of cleaning agent for the time recommended by the detergent's manufacturer.
- ④ Using suitable brushes (never metal brushes or steel wool) cleaning the device thoroughly.
- ⑤ Rinse in running water until all traces of cleaning solution are removed.
- ⑥ Visually inspect for any remaining soil and repeat the steps above if necessary.

2) Cleaning

- ① Immerse the device completely in solution of cleaning agent and activate the ultrasonic bath at the concentration and temperature specified in the detergent manufacturer's instructions.
- ② Follow the same as the clause ④ ~ ⑥ of instructions for pre-cleaning.
- ③ Use a disposable dust-free *absorbent to immediately wipe off the moisture from the instruments

*Absorbent: lint-free and drying material

※ It is recommended to use 0 to 25 us/cm of critical water for the final rinse. (e.g. pre filter, A/C filter, RO filter, 0.22um filter is used or equivalent).

- Sterilization

ARIX Hand / Foot / Wrist / Ankle / Ankle Fibula Hook Plate / Ankle Distal Tibia / Diaphysis / Humerus / Clavicle / Cannulated Screw / Elbow / Small Fragment System must be sterilized before use. Following instructions are recommended;

- Place the devices in the appropriate tray using forceps and/or powder-free gloves to avoid contamination and any other negative effect on the surface of device.
 - Wrap the tray using FDA cleared sterilization wrap (ex. KimGuard KC500).
 - Sterilize in the autoclave validated and maintained in accordance with ISO 17665 and ANSI AAMI ST79. Following parameters are validated in accordance with ISO 17665-1 and recommended for sterilization
- | Cycle | Temperature | Exposure Time | Drying Time |
|------------|-------------|---------------|-------------|
| Gravity | 132°C | 15 min. | 30 min. |
| Pre-vacuum | 132°C | 4 min. | 30 min. |

⚠ Caution: Plates & Screws are disposable (single use only). Do not reuse. Jeil Medical Corporation is not responsible for all problems caused by reuse.

✖ Do not stack trays during sterilization.

[PACKAGING]

All ARIX Hand / Foot / Wrist / Ankle / Ankle Fibula Hook Plate / Ankle Distal Tibia / Diaphysis / Humerus / Clavicle / Cannulated Screw / Elbow / Small Fragment System Plate & Screw are packaged individually. This device is manufactured and sold by Jeil Medical Corporation.

[SYMBOL DESCRIPTIONS]

	Catalogue number
	Batch code
	Do not re-use
	Date of manufacture
	Manufacturer
	Non-Sterile
	Caution
	Consult instruction for use
	Do not use if package is damaged

Rx only United States Federal law restricts medical devices to sale by or on the order of a (licensed healthcare practitioner). (See 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1))

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