JEILMEDICAL

Surgical Instrument

Instruction For Use

Issue Date: 2023.03

Print Date: 2023.03

#837-41000-03

Table of Content

English	03
Español	12
Deutsch	22
Français	33
Italiano	44
Português	55
Polski	66
Αγγλικά	77
Nederlands	88
Čeština	99
Română	108
Magyar	118

[INTENDED PURPOSE]

- Driver Shaft: The driver is intended to enable screw insertion into bone for surgical applications such as orthopaedic surgeries.
- Bender: Used to bend plates to a specific shape.
- Cutter: Used to cut plates or mesh to a specific size.
- Handbody: Handbody is used to drive screws into bone in conjunction with interchangeable driver shafts.
- Holding forceps: The holding forceps is used to hold medical supplies and materials.
- Tray: Reusable sterilization trays for implant and surgical instruments are used for steam sterilization procedures to autoclave surgical instruments in high temperature.
- Depth Gauge: Depth gauge is a measuring device intended for determining the proper length of screws for fastening the ends of a fractured bone.
- Template: It is typically available as a graduated set or an individual plate (e.g., a translucent plastic sheet or a cut metal shape) with patterns/sizes corresponding to the size of the implant to be used to aid the surgeon in final selection of implant size, placement, and/or adaption
- Guide: A manual surgical device designed to simultaneously place, angle, and guide a rotating surgical drill on hard tissue (e.g., bone); some types are also intended for guiding bone-coring instruments (e.g., keel punch).

[MATERIAL]

- Surgical Instrument: Stainless steel etc. (ASTM F899)

[INTENDED USERS]

Surgical instruments are intended to be used in the operating room by dentists and surgeons who are trained, experienced, and familiar with the use of internal bone fixation and techniques. The dentist and surgeon must exercise reasonable judgment when deciding which surgical instrument to use for specific indications.

[PATIENT TARGET GROUP]

This device is used to operate on patients who are applied with the LeForte, Dual top, Leforte Neuro, and ARIX Sternal system. The use of the Leforte system should be considered for its use, indications, taboos, patient bone anatomy, and health conditions.

[CLINICAL BENEFIT]

Bone stabilization and healing are expected to be accelerated due to shortened operation time.

[PERFORMANCE CHARACTERISTICS OF THE DEVICE]

By using this surgical instrument, the operation time is shortened, reducing the patient's surgical delay.

[DESCRIPTION]

Surgical instruments are manufactured and intended for use in surgery using plates and/or screws. The device includes driver shaft, bender, cutter, handbody, holding forceps, tray, and depth gauge. These are non-sterile and prepared in a package and must be sterilized by steam prior to use.

[DIRECTIONS FOR USE]


Surgical instruments offer the operating dentist and surgeon unique features and user-friendly capabilities. A step-by-step explanation of system use is as follows:

1. Before Use


- Make sure that there are no defects on the exterior of the product.
- The product is prepared in a non-sterile state. Please ensure that the product has been properly sterilized before actual use.
- Make sure that it is connected with relevant products.
- The product should be used by a licenced physician (medical specialist).
- Before use, inspect the drill bit for any damage, deformations, and cracks.
- Before surgery, it is necessary to determine which surgical instrument to use based on diagnostic information.

2. How to Use


• Bender


• Cutter


• Driver Shaft


• Tray


• Holding Forceps

• Handbody

2	Insert driver shaft into the protractible end	
3	Drive appropriately for use	

• Depth Gauge

1	Fix the plate and then drill a hole with a drill guide.	
2	Measure the drilling depth by using a depth gauge. Then, insert the appropriate screw into the hook plate hole by using a driver shaft and then fix the screw.	 

• Template

1	Check the shape of the implant to be used in advance.	
---	---	---

• Guide

1	Place the screw or drill bit on the implant to insert it in the correct position.	
---	---	---

[WARNING]

- Ensure that the products have been sterilized properly before use.
- Do not place too much pressure on the products.
- Do not disassemble or modify the product at your own discretion.
- Check the packaging condition for signs of damage.
- Check the assembled condition after assembling instruments and other medical devices.
- The product should be used by a licenced physician (medical specialist).

[CAUTIONS]

- Instructions for use, such as the concentration (0.4 % dilution (Endozime 4 mL/water 1 L)) of the cleaning detergent, should be followed by the cleaning detergent manufacturer's instructions.
- The use of sodium hydroxide is prohibited, and the use of corrosion-causing products or instruments, including rough sponges and metal brushed, should also be avoided.
- It is recommended to use an enzyme cleaner for surgical instruments, and it is not recommended to use a cleaner other than an enzyme cleaner.

[PRECAUTIONS]

- Users should read the instructions before use.
- Users should choose the suitable product model and size before use.
- Please check for any deformations and cracks before use.
- Check the compatibility of instruments before use.
- Do not use products that fall on the floor. If dropped, this product should be sterilized before use.
- Check that the surgical instrument is securely held.
- Pay close attention to patients that may show signs of allergy or hypersensitivity to the chemical components of the material: surgical stainless steel.
- Before each procedure, make sure the pieces are properly seated.

- It is the professional's responsibility to use the Jeil Medical products according to the instructions for use.
- This product should be sterile.
- The non-replacement of the surgical instrument, as recommended by the manufacturer, can generate improper bone heating, compromising the osseointegration success.
- Do not re-grind.
- Do not use solutions with ammonia, acids, or hydrogen peroxide because they can damage the surgical instrument.
- The use of steel brushes is inadvisable. We recommend the use of nylon brushes.
- The improper planning can jeopardize the performance of the implant/prosthesis set causing failures to the system, such as loss or break of the implant, loosening or fracture of the screws.
- Failure to use the drill as recommended by the manufacturer may result in wear, deformation, corrosion, and reduced cutting efficiency.

[CONTRAINDICATIONS]

- Be aware in cases of patients that show signs of allergy or hypersensitivity to the chemical components of the material: Surgical Stainless Steel

[SIDE EFFECTS AND RESIDUAL RISK]

A Residual risk, side effects and adverse effects may occur in all surgical procedures. The dentists and surgeons should operate carefully to prevent problems and take appropriate action for the patient after surgery. Many reactions can occur and the most common are as follows:

Surgical Instrument specific adverse events include but are not limited to:

- Injury to user that may be associated with
 - User damage caused by a sharp edge design
- Infection that may be associated with
 - Use of contaminated product due to sterilization condition error / manufacturing environment contamination / packaging damage / cleaning failure

[Cleaning Procedure]

1. Soaking

- Dilute the cleaning detergent in lukewarm water between (20 to 40°C) to make a suitable cleaning solution, and soak the used instruments for (14 ± 1) minutes. Detachable instruments are disassembled, and then soaked in a cleaning solution. If possible, do not move instruments while cleaning is in operation. For complex structures such as cannulas, grooves, holes, and joints between parts, or in places where the cleaning solution is difficult to reach, use a syringe to spray the cleaning solution several times.

2. Rinsing

- Rinse the instruments under running water between (20 to 40°C) for at least 2 minutes to ensure that no cleaning solution or foreign matter remains.

3. Ultrasonic cleaning

- Prepare a cleaning solution between (20 to 40°C) in the ultrasonic cleaner, and allow the instruments to be completely soaked in the cleaning solution. Clean in an ultrasonic cleaner for 15 minutes.

4. Rinsing

- Rinse the instruments under running water between (20 to 40°C) for at least 2 minutes to ensure that no cleaning solution or foreign matter remains.

5. Ultrasonic rinsing

- Rinse in an ultrasonic cleaner for at least 15 minutes using purified water between (20 to 40°C.)

6. Drying

- Use a disposable dust-free absorbent to immediately wipe off the moisture from the instruments.

[STERILIZATION]

Sterilize in the autoclave validated and follow maintenance procedures in accordance with ISO 17665 and ANSI AAMI ST79. Following parameters are validated in accordance with ISO 17665-1 and recommended for sterilization.

Cycle	Temperature	Exposure Time	Load Characteristics	Drying Time
Gravity	132°C	15 min.	Wrapped*	30 min.
Pre-vacuum	132°C	4 min.	Wrapped*	30 min.

* When handling in accordance with load characteristics, we recommend usage of an FDA-approved wrap to ensure that the device is actually sterile prior to implantation.

[STORAGE AND HANDLING CONDITION]

Careful handling and storage of the product is required. Store in a dry place between 10 to 22°C in temperature and 10 to 30% in humidity.

[DISPOSAL]

This Surgical Instrument is recommended for up to 20 uses, provided that Jeil Medical-recommended conditions of use are followed. Users may use the devices additional times according to their own judgment. A manual to check the condition of the Instrument is provided, and the company is not responsible for additional usage of the instrument at discretion of the user.

[DEVICE Compatibility]

Only compatibility between devices provided by Jeil Medical corporation has been confirmed. Compatibility with other manufacturer's devices has not been tested and Jeil Medical corporation will not possess liability in such instances.

[INCIDENT NOTIFICATION]

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

[SYMBOL DESCRIPTIONS]

SYMBOL	DESCRIPTION	SYMBOL	DESCRIPTION
	Catalogue number		Batch code
	Date of manufacture		Manufacturer
	Non-Sterile		Caution
	Consult instruction for use		Do not use if package is damaged
	Unique device identifier		Medical device
	Temperature limit		Humidity limitation
	Authorized representative in the European Community		EMERGO EUROPE Westervoortsedijk 60, 6827 AT Arnhem, The Netherlands

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